

April 26, 2017

Welcome

Summary

- Cytoscape v3.5.1 is now available for immediate download.
- Addresses 3 issues relative to v3.5 (see Section 10)
 - Help icon on toolbar now loads v3.5 user manual
 - Search bar now finds common words like “a”, “the”, etc.
 - Slider controls now work for small ranges (e.g., 0..1)
- 4 more Apps are available in the App Store as compared to the v3.5 release.
- Ready for use by all segments of the Cytoscape community.
- Continues v3.5’s core feature set, specifically:
 - Upward compatible and can coexist with previous Cytoscape versions
 - Most popular plugins from v2.x versions are available as apps, with more to follow.
- What to do next ...

Release notes: <http://www.cytoscape.org/releasenotes.html>

Dear Cytoscape User

We are proud and pleased to release Cytoscape v3.5.1. This new release contains both new features and improvements to functionality and documentation as described in Sections 7 and 10.

With the v3.5 release, we officially withdraw support for Cytoscape 2.x versions, though installable code and documentation will remain available on the Cytoscape web site.

v3.5.1 retains support for 32 bit JVMs and Windows XP, as they are a significant portion of Cytoscape installations.

v3.5.1 fully supports Java 8, but not Java 7.

This release marks successful tests with the system configurations most common among our user base, including 24 different combinations of operating systems and Java Virtual Machines across Windows, Linux, and Mac OS platforms. Check out your system’s configuration here: <http://tinyurl.com/kelwscg>.

Finally, in v3.5.1 , 4 Apps were added to the App Store (<http://apps.cytoscape.org>), bringing the total to 188. The apps new for v3.5.1 include:

AOP Explorer
cyBrowser
PathInsight
Propagate

We expect to make a v3.6 release in mid Q4 2017. As our plans firm up, we’ll publish them to the Cytoscape community (<http://www.cytoscape.org/roadmap.html>). For more on our roadmap, see Section 9.

If you are a new or casual Cytoscape user, **welcome aboard!** Feel free to either read on or jump ahead to *What to Do Next* for a quick start.

For Cytoscape v2.x Users

As of the Cytoscape v3.5 release on March 27, 2017, the Cytoscape releases at the 2.x level are no longer supported. This reflects both the great success Cytoscape v3 has had since its initial release on February 4, 2013 and the rapidly diminishing user base for Cytoscape 2.x.

While all version 3 releases are upwardly compatible with previous Cytoscape versions (i.e., they can load, analyze, and display networks and attributes stored using previous versions, including v2.x), they can coexist with previous versions on the same workstation. It is not necessary to uninstall previous versions before installing a newer release, and you can choose to use both on the same workstation.

Version 3 has matured and achieved a reputation for speed, reliability, and support that makes it a *superior choice* (compared to v2.x) for most uses. If you are using v2.x, please switch to v3.5 and enjoy the new v3 features and the Apps available in the App Store (<http://apps.cytoscape.org>).

A table describing Cytoscape support is presented in Section 8.

We are very excited to make this release, and hope you find it an indispensable tool in your research. To get you started, we have provided additional introductory material to set you on the right foot, and we have also provided background material for general interest.

Please let us know how we're doing and how we can better advance your research!

The Cytoscape Team

What to Do Next

We hope to help you use Cytoscape in the most productive way. Your next steps depend on how experienced you are with Cytoscape and what you would like to do. In the table below, please find the category that best describes you, and proceed accordingly.

Type of Cytoscape User	Where to Read Next
New	For New Users (Section 1)
Casual	For Casual Users (Section 2)
Power	For Power Users (Section 3)
App Developers	For App Developers (Section 4)
General Interest	Try getting oriented by reading: General Background (Section 5) Cytoscape App Store (Section 6) Then try browsing: New Features (Section 7) Support Offering (Section 8) Future Plans (Section 9) To get your feet wet, try: For New Users (Section 1)

1. For New Cytoscape Users

What is Cytoscape?

Cytoscape is software that can help you visualize complex networks (as *graphs*) representing relationships between genes, proteins, or other entities. Cytoscape tracks and displays these entities (as graph *nodes*) and relationships between them (as graph *edges*). Cytoscape, per se, does not assign meaning to either nodes or edges – you can assign your own meanings to them, and you can add your own node or edge attributes (e.g., name, measurements, states, etc.) that customize your graph to represent semantics important to you. While Cytoscape provides simple analytics and visualizations, others are available as apps from the Cytoscape App store (see Section 6).

What Do I Need Before I Can Use Cytoscape?

At the very least, you must provide Cytoscape with a list of nodes. You can also provide a list of edges and attributes. You can enter your graph manually or import it from your own files – Cytoscape supports a wide variety of file formats, including text and Excel. Cytoscape can also load pre-configured networks (e.g., from BioGRID) so you can leverage them without having to enter them yourself.

Can Cytoscape Help Me Get Started?

Yes! The Cytoscape team has provided many resources tailored for all kinds of users (see Section 8). To start, there is an Installation Guide and numerous tutorials. Once you load your graph, you may find additional useful analytics and visualizations at the App Store (see Section 6).

Which Cytoscape Version Should I Use?

Generally, you should use v3.5.1 (the newest version). v2.x has been in the field for several years and is now out of support. Cytoscape v3.5.1 has many improvements over previous Cytoscape v3.x releases.

Note that if you have already installed a previous Cytoscape version, you can install v3.5.1 without uninstalling your previous version.

How to Download and Install Cytoscape

- Surf to <http://www.cytoscape.org>, click on the Download Cytoscape button, and follow the prompts
- ***Be sure to consult with the Release Notes (viewable via the Release Notes button) for issues and solutions pertaining to your installation***
- Surf to <http://manual.cytoscape.org> and follow the instructions in the *Launching Cytoscape* section.

How to Troubleshoot Cytoscape Installation?

Normally, the main requirement for installing Cytoscape is a Java Virtual Machine (JVM) appropriate for your workstation. For Windows users, the JVM commonly installed is a 32 bit version, even if the Windows itself is 64 bits. ***For Cytoscape use, the 64 bit version is highly recommended instead.*** For all workstations, Java 8 is required.

Regardless of your workstation operating system, you can learn whether your environment is sufficient to execute Cytoscape by following the Troubleshooting instructions at <http://www.cytoscape.org/troubleshooting.html>.

2. For Casual Cytoscape Users

If you're confident of your Cytoscape skills, please read the Power Users information (Section 3).

If you're Cytoscape skills are rusty, please continue with the New Users information (Section 1).

3. For Cytoscape Power Users

Is v3.x for Me?

The answer is most likely yes! Cytoscape 2.x is out of support.

Most modern v2.x plugins have been converted to v3.x apps, and new v3.x apps are more plentiful than v2.x plugins (-- *app* is the v3.x name for the installable features known as a *plugin* in v2.x). To take stock, you can browse the App Store directly (Section 6).

Note that v2.x plugins are not compatible with v3.x apps. If you rely on v2.x plugins that are not yet available for v3.x, you may need to remain with v2.x instead of upgrading to v3.x, though v2.x is no longer supported.

Can I use v2.x and v3.x at the Same Time?

You can experiment with v3.x by installing it on your workstation – you can use both v2.x and v3.x on the same workstation.

How to Prepare for Cytoscape

1. If you have already installed a v3.5.1 RC version, *please uninstall it before proceeding*. Previous Cytoscape versions can coexist with v3.5.1, and do not need to be uninstalled.
2. Consult the App Store (Section 6) to verify that the apps you need are available. Some apps have become part of the Cytoscape core, and you won't need to download them. Other apps are available in the Cytoscape App store.
3. Browse the Future Plans table (Section 9) and Support Offering table (Section 8) to see where the Cytoscape community is going.

How to Download and Install Cytoscape

- Surf to <http://www.cytoscape.org>, click on the Download Cytoscape button, and follow the prompts
- ***Be sure to consult with the Release Notes (viewable via the Release Notes button) for issues and solutions pertaining to your installation***
- Surf to <http://manual.cytoscape.org> and follow the instructions in the *Launching Cytoscape* section.

How to Troubleshoot Cytoscape Installation?

Cytoscape v3.5.1 requires Java 8, preferably with 64 bit support. You can learn whether your environment is sufficient to execute Cytoscape by following the Troubleshooting instructions at <http://www.cytoscape.org/troubleshooting.html>.

4. For Cytoscape App Developers

Cytoscape manages graph storage and display, and contains a number of common analytic and visualization functions. App Developers may write *apps* (formerly known as *plugins*) to perform problem-specific analytics, graph layouts, or other visualizations.

Before starting app development, you should be familiar with overall Cytoscape usage and functionality – see the Power Users section (Section 3) to start working with Cytoscape if you have not already done so.

Note that v2.x plugins do not work with v3.x. In order to get the benefits of improved stability, the Cytoscape team found it necessary to change the relationship between plugins/apps and the Cytoscape core. Many plugins have been ported to v3.x (see Section 6). The porting exercise is substantial, and once a plugin is converted to an app, it will continue to work at least until v4.0 is released – there is no v4.0 planned at this time. For more on this topic, see:

http://wiki.cytoscape.org/Cytoscape_3/CoreDevelopment/Motivation

Should I Convert My v2.x Plugin to a v3.x App?

Yes! Cytoscape v2.x is now out of support.

How Do I Get Started With App Development?

To get started with app development:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/Cytoscape_App_Ladder

To convert a v2.x plugin to a v3.x app, refer to:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/PluginPortingGuide

or

<http://f1000research.com/articles/3-150/v1>

For a broader view of the Cytoscape system (including the Core), refer to:

http://wiki.cytoscape.org/Cytoscape_3/CoreDevelopment

Are there App Samples?

Yes! See:

http://opentutorials.cgl.ucsf.edu/index.php/Tutorial:Creating_a_Simple_Cytoscape_3_App

Also note the App Cookbook:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/Cytoscape_3_App_Cookbook

5. General Background

The overall mission of Cytoscape is to be a freely available worldwide asset supporting network analysis and visualization for systems biology science.

The major focus of v3.x is the modularization and rationalization of code to solve stability issues in v2.x encountered as multiple developers pursued multiple agendas. Under v2.x, internal programmatic interfaces evolved from one release to the next, leading to the failure of working plugins over time and negative interactions between otherwise working plugins. Ultimately, this resulted in loss of programmer and user productivity, and undermined community confidence in Cytoscape.

v3.x addresses these issues by adopting modular coding practices promoted by the OSGi architectural framework¹. This enables both the Cytoscape core and externally developed apps (formerly called *plugins*) to evolve independently without compromising unrelated functionality. At the logical level, Cytoscape leverages OSGi precepts to produce v3.x APIs having cleaner and clearer demarcations between functional areas. At the deployment level, OSGi enables on-the-fly substitution of one processing element for another (e.g., apps) in order to tailor Cytoscape to meet user requirements at runtime without reinstalling or reconfiguring Cytoscape.

Creating v3.x occupied an international team of engineers for over a year, and represents a strong investment toward reducing future development and support costs, and increasing reliability and evolvability. With the completion of v3.x, we expected to leverage v3.x as a platform to satisfy the evolving needs of multiple stakeholder groups, and as a platform enabling research on leading edge analysis and visualization techniques. v3.x is the successor to v2.x, with 2.x now out of support.

During the creation of v3.x, a number of features were added, as listed in the New Feature Table below.

¹ www.osgi.org – also used as the basic framework for Eclipse and numerous commercial products

6. The Cytoscape App Store

The Cytoscape App Store targets broad demographics through specific features:

Cytoscape User	Benefit
New or Casual	Find Cytoscape functionality easy via familiar paradigms of search and tag navigation
Power or Frequent	Find and install apps easy and provide crowdsourced rating and reviews
App Developers	Have a dedicated presence for their work, including a direct line to their users and critical statistics to justify their development on the Cytoscape platform.

Furthermore, the App Store supports contests and social incentives that encourage critical technical and community development, including porting, analytics, rating/reviewing, and so on.

Note that the App Store supports only Cytoscape v3.x and no longer accepts plugins compatible with earlier Cytoscape versions. Likewise, plugins that work with earlier versions are not compatible with v3.x – the App Store (<http://apps.cytoscape.org>) contains a list of apps available for v3.x.

To access the App Store from Cytoscape, access Cytoscape's new Apps menu and choose the App Manager menu item.

7. New Feature Table

Feature	Beneficiary	Manual Section
ID Mapper	All users	8.4. Mapping Identifiers
Export Dialog Improvements	All users	22. Export Your Data
CoSE Layout (for compound nodes)	Power users	12.3. Automatic Layout Algorithms
New API (CyJob) for Long Running Tasks	App developers	
Drag/Drop for Network and Table Files	All users	5.1. Import Fixed-Format Network Files 8.1. Import Data Table Files
Consolidated Node Layout Tools	Power users	12.8. Node Layout Tools

8. Support Offering Table

Feature	Purpose	Beneficiary	URL
Cytoscape Google Groups <ul style="list-style-type: none">HelpdeskDiscuss	Community response on best effort basis	All users	mailto:cytoscape-helpdesk@googlegroups.com mailto:cytoscape-app-dev@googlegroups.com
User Manual	Comprehensive manual	All users	http://manual.cytoscape.org
Installation Guide	Cytoscape installation	New/Casual users	Section 2 of http://manual.cytoscape.org
Tutorials	Streamlined productivity	All users	http://opentutorials.cgl.ucsf.edu/index.php/Portal:Cytoscape3
Presentation Slides on Web Site	Enable casual learning	All users	https://speakerdeck.com/cytoscape
Code Snippets	Commonly used App development code	App Developers	http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/Cytoscape_3_App_Cookbook

9. Future Plans Table

Feature	Purpose	Beneficiary	Timeframe
Distributed computing	Link to NDEx and other external storage and compute services	All users	Ongoing
Additional Visualizations			Ongoing
Additional Analytics			Ongoing
Wallpaper Collection	Publicize good results associated with Cytoscape	All users, potential users, general public	TBD
Cytoscape Clinics New/Casual User Power User App Developer	Chat rooms give instant response based on community availability	All users	Ongoing

10. Issues Addressed

The following issues were addressed in v3.5.1. Details for each issue can be found at Cytoscape's RedMine issue repository:

<http://code.cytoscape.org/redmine>

Issue	Description
3808	Remove Lucene stop words from search
3807	Help button (top left of worksheet) goes to outdated link
-	Slider does not work for narrow ranges